



**INSTRUMENT STRINGERS**


**SteriBest® Instrument Stringers** are manufactured from medical grade, high quality 300 series stainless steel in a wide variety of beneficial lengths, widths, and styles to meet the needs of virtually every medical facility.

- Manufactured with high quality 300 series stainless steel.
- Many stringer styles available.
- Most styles are available in two widths.
- Private label programs and custom products available.



**CENTER LATCH STRINGER**  
3.0 in. width

STCL-6	6"	STCL-14	14"
STCL-8	8"	STCL-16	16"
STCL-10	10"	STCL-18	18"
STCL-12	12"		



**CENTER LATCH STRINGER**  
5.75 in. width

STCL55-6	6"	STCL55-12	12"
STCL55-8	8"	STCL55-14	14"
STCL55-10	10"	STCL55-16	16"



**PIN LOCK STRINGER**  
2.75 in. width

ST-4	4"	ST-12	12"
ST-6	6"	ST-14	14"
ST-8	8"	ST-16	16"
ST-10	10"		



**PIN LOCK STRINGER**  
3.0 in. width

ST-18	18"
ST-20	20"



**PIN LOCK STRINGER**  
5.5 in. width

ST55-6	6"	ST55-12	12"
ST55-8	8"	ST55-14	14"
ST55-10	10"	ST55-16	16"



The use of instrument stringers is integral to the preparation of instruments for sterilization. Stringers hold ring-handled instruments open during the cleaning process, to ensure the maximum surface area of the instruments are exposed to the cleaning solution. A closed instrument may conceal bioburden which, if left intact, would require additional cleaning methods. The open position of the instruments on the stringers also ensures better exposure to heated air drying, for quicker dry times.




**OPEN END STRINGER**  
2.75 in. width

U-OE-ST4	4"	U-OE-ST12	12"
U-OE-ST6	6"	U-OE-ST14	14"
U-OE-ST8	8"	U-OE-ST16	16"
U-OE-ST10	10"	U-OE-ST18	18"



**OPEN END STRINGER**  
5.5 in. width

U-OE-ST5-6	6"	U-OE-ST5-12	12"
U-OE-ST5-8	8"	U-OE-ST5-14	14"
U-OE-ST5-10	10"	U-OE-ST5-16	16"



**SPRING STRINGER**  
2.75 in. width

SPS-6	6"	SPS-12	12"
SPS-8	8"	SPS-14	14"
SPS-10	10"	SPS-16	16"



**T-BAR STRINGER**

STB-6	6" length
STB-8	8" length
STB-12	12" length



**WEINSTEIN STRINGER**  
3.0 in. width

WSR-6	6"	WSR-12	12"
WSR-8	8"	WSR-14	14"
WSR-10	10"	WSR-16	16"



**INSTRUMENT STAND**

STR-STAND3.5

The instrument stand holds the instruments upright and open. The stand can be secured to a basket or remain free. Comes with four 3.5mm screws for securing to a basket.


Instrument stringer width measured outer edge to outer edge.

## INSTRUMENT STRINGERS



**PIN AND SOCKET STRINGER**  
3.0 in. width

STPS-6	6" length
STPS-8	8" length
STPS-10	10" length
STPS-12	12" length



**INSTRUMENT STRINGER HOLDING BOARD**  
7.75" X 10.5"

SHB-B1  
Holds stringers while instruments are being placed on them. (*Stringer not included.*)



**BUNT FORCEPS HOLDER**

BUNT-FH



**BAND CLAMPS**

H-40G	40 mm
H-50G	50 mm
H-60G	60 mm
H-70G	70 mm
H-80G	80 mm

Made of stainless steel and coated with green silicone.  
Allows instrument jaws and joints to remain open during cleaning.



**INSTRUMENT SPREADER**


H-00

Made of stainless steel.  
Keeps ring-handled instrument jaws open during cleaning.

# ADJUSTABLE INSTRUMENT STRINGERS



Adjustable Stringers are designed for use in any width or length basket or tray allowing for stringing a variety of open instrument sizes. Stringers adjust in width and are available in three lengths.



**ADJUSTABLE STRINGER**  
3.5 in. - 5.5 in. width

**ADJUSTABLE STRINGER**  
3.5 in. - 5.5 in. width

**ADJUSTABLE STRINGER**  
3.5 in. - 5.5 in. width

AJS8      8" length

AJS13      13" length

AJS17      17" length