




Weck® Hemoclip® Ligation Clip System

Hemoclips offer veterinarians several advantages over sutures:

- Significant reduction in surgical time.
- Improved access in confined and hard-to-reach surgical areas.
- Improved surgical site visibility.
- Reduced patient table time resulting in lower anesthesia and overall operating room costs.



It has been documented that Hemoclips will save veterinarians substantial time and money, especially when used in routine surgical procedures.

The ligation of blood vessels proceeds faster than the time to clamp and ligate vessels with conventional suture material. This results in a significantly shorter surgical time for the patient and the veterinarian. In one study, the use of Hemoclips in ovariohysterectomies reduced the surgical time by approximately 24% in dogs and 33% in cats (Table 1).

TABLE 1*

	Average Surgical Time for Ovariohysterectomy (min.)		Average Reduction In Surgical Time
	Without Hemostatic Clips	With Hemostatic Clips	
DOGS	17.4 minutes	13.3 minutes	24%
CATS	13.2 minutes	8.9 minutes	33%

Veterinarians report Hemoclips are of great value for vessel ligation in the following types of surgical procedures:

- Feline and Canine Ovariohysterectomy
- Feline, Canine and Equine Castration
- Surgical Removal of Cutaneous and Subcutaneous Neoplasms
- Mastectomy
- Splenectomy, Nephrectomy, Adrenalectomy, Thyroidectomy
- Limb Amputation
- Tail Amputation, Toe Amputation
- Ocular Enucleation
- Intestinal Resection
- Patent Duct Arteriosus (PDA)
- Canine Vasectomy
- Liver Biopsy
- Persistent Right Aortic Arch
- Avian and Exotic Animal Surgical Procedures
- Equine Colic Surgery
- Plus many, many more...

See reverse side for the actual sizes of the sterile, titanium Hemoclip clips and current pricing.

* Zymet, C.L.: A Rapid Modified Approach to Feline and Canine Ovariohysterectomies Using Hemostatic Clips. *Journal of American Animal Hospital Assoc.* 10: 566-569, 1974.

Weck® Hemoclip® Ligation Clip System

2022 Veterinary Price List

LIST #	HEMOCLIP SIZE	CLIPS/ CARTRIDGE	CARTRIDGES/ BOX	CLIPS/ BOX	VETERINARY PRICE PER BOX	VETERINARY PRICE PER CLIP
--------	---------------	---------------------	--------------------	---------------	-----------------------------	------------------------------

Traditional Hemoclip System – Titanium Clips

523735	Small (yellow)	25	12	300	\$266.90	\$0.89
523700	Medium (blue)	25	10	250	\$262.95	\$1.05
523760	Medium/Large (green)	10	20	200	\$268.50	\$1.34
523770	Large (orange)	15	10	150	\$270.35	\$1.80

523835	Small, Ten Pak (yellow)	10	24	240	\$234.95	\$0.98
523800	Medium, Ten Pak (blue)	10	20	200	\$231.10	\$1.16
523860	Medium/Large, Ten Pak (green)	10	16	160	\$236.90	\$1.48
523870	Large, Ten Pak (orange)	10	12	120	\$237.60	\$1.98


LIST #	SAMUELS HEMOCLIP APPLIERS	LENGTH	TYPE	PRICE
523140	Samuels Applier, Small	6.0"		\$172.50
523109	Samuels Applier, Medium	6.0"		\$172.50
523105	Samuels Applier, Medium	8.0"	Right Angled	\$172.50
523165	Samuels Applier, Medium/Large	8.0"		\$172.50
523167	Samuels Applier, Medium/Large	8.0"	Right Angled	\$172.50
523175	Samuels Applier, Large	8.0"		\$172.50
523177	Samuels Applier, Large	8.0"	Right Angled	\$172.50

LIST #	CARTRIDGE BASE TO ASSIST LOADING	LENGTH	PRICE
523130	Hemoclip Ligation Cartridge Base	7.0"	\$180.00

DIMENSIONS OF WECK HEMOCLIPS

HEMOCLIP DESCRIPTION	WIDTH		HEIGHT		CLOSED CLIP LENGTH	
	INCHES	MILLIMETERS	INCHES	MILLIMETERS	INCHES	MILLIMETERS
Small Hemoclips (yellow)	0.085	2.16	0.105	2.67	0.140	3.55
Medium Hemoclips (blue)	0.142	3.59	0.177	4.49	0.228	5.79
Med-Large Hemoclips (green)	0.210	5.33	0.296	7.51	0.365	9.27
Large Hemoclips (orange)	0.251	6.37	0.415	10.53	0.500	12.71

NOTE: Width and height are taken from **INSIDE** the clip. For informational use only, all dimensions are **APPROXIMATE**.



To order, contact your authorized RICA Veterinary Distributor for current pricing, call RICA's office in Schiller Park, Illinois at 1-800-889-3218 or visit www.ricasurgical.com.



RICA Surgical Products, Inc.
 9207 Ivanhoe Street
 Schiller Park, IL 60176
 (800) 889-3218



PRINTED IN USA
 FORM # RSP-WHCLCS-0122