

# STERIBEST<sup>®</sup>.com

BETTER HEALTH THROUGH BETTER CIRCULATION

## Micro Instrument Trays from RICA Surgical Products



**SteriBest<sup>®</sup> Micro Instrument Trays** offer advanced design features which facilitate flow and penetration of Steam and ETO as well as promote faster, more complete drying, as demonstrated in independent lab studies.

**Pin Mats** and **Silicone Bar Holders** are designed to minimize instrument movement, protect instruments, and maximize sterilant flow and contact.

**Stackable** for storage.

**Made in USA.**

### Independent Test Lab Validations:

#### • DRY TIME STUDY


**Dry Time Validation** at 25 minutes, tray is completely dry both above and below the mats, helps you establish your parameters to ensure sterile instruments are moisture free.

#### • STERILIZATION VALIDATION STUDY

**Sterilization Validation** shows 100% kill both above and below the mats at 2 minutes to give you confidence that you are delivering sterile instruments to the procedure.



Refer to our website for full studies: [Steribest.com](http://Steribest.com)

# Steribest Instrument Trays



## Item CP1038D2

Base, Lid, Mat, Insert Tray, Insert Mat  
10.0 in. x 6.0 in. x 1.5 in.  
254mm x 152mm x 38mm



## Item CP1038D1

Base, Lid, Mat  
10.0 in. x 6.0 in. x 1.5 in.  
254mm x 152mm x 38mm



## Item CP1038D1-P-1S/L

Base, Lid  
10.0 in. x 6.0 in. x 1.5 in.  
2 Bars with 1-Slot  
Phaco Holder  
16mm & 8mm


#CP1038D1-P-1S/L



## Item CP1038D1-P-3S/L

Base, Lid  
10.0 in. x 6.0 in. x 1.5 in.  
2 Bars with 3-Slot  
Phaco Holder  
One bar has three 8mm slots;  
Other bar has three slots:  
15.5mm, 7.5mm, and 12.5mm


#CP1038D1-P-3S/L



## Item CP1038S Base, Lid, Mat

10.0 in. x 6.0 in. x .75 in.  
254mm x 152mm x 19mm

Item #CP1038S



Item #CP471D1

**Item CP471D1**

Base, Lid, Mat  
 7.5 in. x 4.0 in. x 1.5 in.  
 190.5mm x 102mm x 38mm



#CP471D1-P-1S/L

**Item CP471D1-P-1S/L**


Base, Lid  
 7.5 in. x 4.0 in. x 1.5 in.  
 2 Bars with 1-Slot  
 Phaco Holder  
 16mm and 8mm



#CP471D1-P-3S/L

**Item CP471D1-P-3S/L**

Base, Lid  
 7.5 in. x 4.0 in. x 1.5 in.  
 2 Bars with 3-Slot Phaco Holder  
 One bar has three 7.5mm slots;  
 Other bar has three slots:  
 15.5mm, 7.5mm and 12.5mm



Item #CP471S

**Item CP471S**

Base, Lid, Mat  
 7.5 in. x 4.0 in.  
 x .75 in.  
 190.5mm x 102mm  
 x 19mm



#CP471S-DK-2

**Item CP471S-DK-2**

Base, Lid  
 7.5 in. x 4.0 in. x .75 in.  
 2 Bars with 2-Slot  
 Knife Holder  
 5.5mm



#CP471S-DK-3

**Item CP471S-DK-3**

Base, Lid  
 7.5 in. x 4.0 in. x .75 in.  
 2 Bars with 3-Slot  
 Knife Holder  
 5.5mm



#CP471S-DK-4

**Item CP471S-DK-4**

Base, Lid  
 7.5 in. x 4.0 in. x .75 in.  
 2 Bars with 4-Slot  
 Knife Holder  
 5.5mm



#CP471S-DK-5

**Item CP471S-DK-5**

Base, Lid  
 7.5 in. x 4.0 in. x .75 in.  
 2 Bars with 5-Slot Knife Holder  
 5.5mm

# Instrument Trays (continued)



Item #CP614

## Item CP614

Base, Lid, Mat  
6.0 in. x 2.5 in. x .75 in.  
152mm x 64mm x 19mm



## Item CP614-DK-1

Base, Lid  
6.0 in. x 2.5 in. x .75 in.  
2 Bars with 1-Slot  
Knife Holder  
5.5mm

#CP614-DK-1



## Item CP614-DK-2

Base, Lid  
6.0 in. x 2.5 in. x .75 in.  
2 Bars with 2-Slot  
Knife Holder  
5.5mm


#CP614-DK-2



## Item CP614-DK-3

Base, Lid  
6.0 in. x 2.5 in. x .75 in.  
2 Bars with 3-Slot  
Knife Holder  
5.5mm

#CP614-DK-3



Item #CP634

## Item CP634


Base, Lid, Mat  
6.5 in. x 4.0 in. x .75 in.  
165mm x 102mm x 19mm



#CP634-DK-5

## Item CP634-DK-5

Base, Lid  
6.5 in. x 4.0 in. x .75 in.  
2 Bars with 5-Slot Knife Holder, 5.5mm



Item #CP778

## Item CP778

Base, Lid, Mat  
7.5 in. x 2.5 in. x .75 in.  
190.5mm x 63.5mm x 19mm



## Item CP778-DK-1

Base, Lid  
7.5 in. x 2.5 in.  
x .75 in.  
2 Bars with  
1-Slot  
Knife Holder  
5.5mm

#CP778-DK-1



## Item CP778-DK-2

Base, Lid  
7.5 in. x 2.5 in.  
x .75 in.  
2 Bars with  
2-Slot Knife  
Holder  
5.5mm

#CP778-DK-2



## Item CP778-DK-3

Base, Lid  
7.5 in. x 2.5 in.  
x .75 in.  
2 Bars with  
3-Slot  
Knife Holder  
5.5mm

#CP778-DK-3